TECHNICAL BULLETIN

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Date of the document issue: the 20th of April 2018

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Gel based on modified hyaluronic acid hydrous sterile

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2.0 ml of HYALREPAIR[®]-02 Chondroreparan[®] in a glass syringe with Luer-Lock tip. The content is sterile. Syringe with gel, 2 adhesive labels for physician attached to a syringe are enclosed in a sealed blister. One blister is packaged with instruction for use in a cardboard box.

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Individual package – cardboard box, size 30 x 77 x 197 mm.

Group package – 42 individual packages in a corrugated cardboard box, size 309 x 302 x 257 mm.

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HYALREPAIR[®]-02 Chondroreparan[®] contains a unique form of high-purity modified hyaluronic acid obtained by bacterial fermentation with addition of hyaluronic acid copolymer with magnesium ascorbylphosphate, L-proline, L-lysine and glycine. It creates optimum environment for normal functioning of fibroblasts which are responsible for synthesis and secretion of intercellular matrix proteins (collagen, elastin) and mucopolysaccharides.

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HYALREPAIR[®]-02 Chondroreparan[®] is used in sports medicine, rheumatology, traumatology, orthopedics and maxillofacial surgery. It may be administered into any joints.

Indications:

1. As a part of the combined therapy:

– treatment of I-III stage osteoarthritis of knee, coxofemoral and other synovial joints of degenerative-dystrophic or traumatic origin;

– necessity to improve viscoelastic and protective properties of the synovial fluid;

– degenerative-dystrophic conditions of the ligamentous apparatus, connective tissue structures and cartilages;

– arthritis in remission.

2. For prevention of degenerative-dystrophic diseases of the musculoskeletal system.

3. During the rehabilitation period after injuries and surgical interventions.

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The amount of the injected gel and the type of a needle are determined by a physician on a per-patient basis depending on the anatomic and physiological characteristics and the condition of the joint.

In case of effusion removal of the excessive synovial fluid from the joint before the injection of the gel is recommended.

Depending on the size of the joint, it is necessary to inject into the joint up to 2 ml of HYALREPAIR[®]-02 Chondroreparant[®] every 14 days for 6-8 weeks. Refresher course is conducted in 6-12 months. It is possible to treat several joints simultaneously.

Treatment with HYALREPAIR[®]-02 Chondroreparant[®] produces the effect only on the joint into which the gel is injected, HYALREPAIR[®]-02 Chondroreparant[®] does not produce a systemic effect.

Important! It is vital to comply with the injection guide to ensure due result of the procedure. Injection of HYALREPAIR[®]-02 Chondroreparant[®] may only be performed by specialists who have been duly authorized there to according to the local legislation and properly trained.

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To store in the place protected from the sun light at a temperature of +5 to +15 °C. Do not freeze. Do not heat.

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Production is manufactured in accordance with ТУ 9398-004-58568834-2009.

Quality Management System of the Company-manufacturer is certified according to the requirements of EN ISO 13485:2012.

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